

# Experience your perfect combination

In only 4 steps. Choose your perfect configuration.

1



S-FRAME



Professional L Headband

## MOUNTING OPTION

In a lightweight professional design for fatigue-free work

2



HR 2.5x/340mm



HR 2.5x/420mm



HR 2.5x/520mm



HRP 3.5x/420mm



HRP 4x/340mm



HRP 6x/340mm

## MAGNIFICATION AND WORKING DISTANCE

Choose the smallest magnification that fulfils your requirements. The lower the magnification, the easier it is to use the loupes due to the larger field of view and better depth of focus. The working distance depends on personal preference.

3



LED NOW IN HERE QUALITY.

## ILLUMINATION

Choose between the LoupeLight2 and the ML4 LED.



LoupeLight2  
for S-FRAME

Compact and lightweight.  
Coaxial illumination.



ML4 LED HeadLight  
for Professional L Headband

Coaxial illumination.  
Flexible power source options.

4

## POWER SOURCE



mPack mini  
Portable  
power pack



mPack  
Portable  
power pack



mPack UNPLUGGED  
Wireless.  
Headband-mounted

For more  
information,  
please  
contact us!